

Blood Donation/ Collection Chair

Electric Version - Série IV



The Blood Donation/Collection Actualway® Série IV Chair is intended to be used in laboratories, clinics, blood services and can also be used for small surgeries, gynaecological examination and recovery

FEATURES AND ADVANTAGES


- Easy to use and maintain
- Excellent delivery date and efficient after sales service
- 3 motors: back rest, height and leg rest regulation
- Adjustable head cushion
- Handset
- Adjustable, rotating and removable arm rests which can be used as leg holders
- Granted stability
- Tested and approved loaded capacity: 200Kg + Accessories
- Rotating rear castors, fixed front castors
- Individual locking system on rear castors
- 10cm diameter castors
- High quality electrical components and accessories
- Comfortable upholstery with appropriate density
- Class M2 fireproof and waterproof French coated fabric, efficient against bacteria (MRSA), easy to clean with water and soap, with medical properties
- Wide range of available colours

DIMENSIONS

- Min. height: 50cm
- Max. height: 70cm
- Overall width: 88cm
- Overall length: 175cm
- Seat rest width: 54cm
- Seat rest: 48x54x10cm
- Back rest: 76x54x10cm
- Leg rest: 45x54x10cm
- Arm rest: 48x15x4,5cm
- Base metal frame: 95x62x25cm

LOGISTIC DATA

- Net weight: 75Kg
- Pallet: 0,80x1,20x1,60 m
- Gross weight: 100Kg




LARGE RANGE OF OPTIONS

- Button to get to the horizontal position, max. -4°
- Stainless steel IV serum pole with 4 hooks
- Stainless steel paper roll holder
- Stainless steel tray with support
- Seamless upholstery
- Class M1 fireproof and waterproof coated fabric (Urban, Sirius...)
- Back handle to move the chair easily
- Special arm rest support regulation and fixation
- 110V control box
- Battery

TECHNICAL SPECIFICATIONS

- Power supply: 230V / 50-60Hz
- Protection Classification: IP44
- Handset Protection Classification: IP66

Distributor:



Producer:
VIA ACTUAL-COMÉRCIO INTERNACIONAL, LDA
CE Castelo da Maia
Rua Manuel Assunção Falcão, 481- Fracção H
4475-041 Santa Maria de Avioso - Maia
PORTUGAL

Tif. / Phone: +351 229 829 663/4
Fax: +351 229 829 665
E-mail: comercial@actualway.pt
www.actualway.pt