

# Blood Donation/ Collection Chair

## Manual Version Gas spring - Série VI



*The Manual Blood donation/collection Actualway® Série VI Chair is intended to be used in laboratories, clinics and blood services*

### FEATURES AND ADVANTAGES

- Easy to use and maintain
- Excellent delivery date and efficient after sales service
- Multi positions
- 3 sections
- Synchronized movements of the back rest and leg rest
- Gas-spring operated
- Adjustable, removable and rotating arm rests
- 4 castors
- Individual locking system on the rear castors
- Rotating rear castors, fixed front castors
- 10cm diameter castors
- Adjustable head cushion
- Granted stability
- Tested and approved loaded capacity: 175Kg + Accessories
- Comfortable upholstery with appropriate density
- Class M2 fireproof, waterproof French coated fabric, efficient against bacteria (MRSA), easy to clean with water and soap, with medical properties
- Wide range of available colours


### DIMENSIONS

- Width: 54cm
- Overall length: 171cm
- Overall width (including arm rest): 90cm
- Seat rest: 50x54x10cm
- Back rest: 76x54x10cm
- Leg rest: 30x54x10cm
- Arm rest: 40x16x5cm
- Base metal frame: 95x62x25cm

### OPTIONS


- Stainless steel IV serum pole with 4 hooks
- Stainless steel paper roll holder
- Stainless steel tray with support
- Seamless upholstery
- Class M1 fireproof and waterproof coated fabric (Urban, Sirius...)
- Back handle to move the chair easily
- Arm rest support regulation and fixation mechanism

Distributor:



### LOGISTIC DATA

- Net weight: 65Kg
- Pallet: 0,80x1,20x1,60m
- Gross weight: 85Kg



Producer:  
VIA ACTUAL-COMÉRCIO INTERNACIONAL, LDA  
CE Castelo da Maia  
Rua Manuel Assunção Falcão, 481- Fracção H  
4475-041 Santa Maria de Avioso - Maia  
PORTUGAL

Tlf. / Phone: +351 229 829 663/4  
Fax: +351 229 829 665  
E-mail: comercial@actualway.pt  
www.actualway.pt