

# Hemocompact Treatment Chair Série V

The best way to take care



The Hemocompact Actualway® Chair Série V was developed to be functional and versatile and is intended to be used in Dialysis, Small Surgeries, Recovery, Observation, Sampling and Blood Transfusions

## FEATURES

- Easy to use and maintain
- Special dimensions for small treatment rooms
- Excellent delivery date and efficient after sales service
- 3 motors (back rest, height and leg rest adjustment)
- Adjustable head cushion
- Handset including the Trendelenburg operating key
- Arm rest inclination adjustment
- Granted stability
- Tested and approved loaded capacity: 200Kg + Accessories
- Rotating rear castors
- Individual blocking system
- 10 cm diameter castors
- Included battery
- High quality electrical components and accessories
- M2 Fireproof and waterproof coated fabric, efficient against bacteria (MRSA)
- Wide range of available colours

## LARGE RANGE OF OPTIONS

- Weight scale M and Class III
- Built-in emergency crash button very easy to activate to get the Trendelenburg position (-12°) in 12 sec.
- Optional big leg rest
- MCL control box: 110V
- IV serum pole
- Back Handle
- Paper roll holder
- Foldable and adjustable, 5Kg capacity, rectangular ABS tray
- Upholstery colour range class M1 (Urban, Sirius...)
- Seamless

## TECHNICAL SPECIFICATIONS

- Power supply: 230V /50-60 Hz
- Protection Classification: IP44
- Handset and MCL protection category: IP66

## LOGISTIC DATA

- Net wet: 75Kg
- Gross weight: 100Kg
- Pallet: 0,80x1,20x1,60m

## DIMENSIONS

- Min./Max. height: 55-75 cm
- Dimensions on the standard position: 62x110x117cm
- Overall length: 175cm or 187cm (with optional big leg rest)
- Width: 54cm
- Total width: 88cm
- Frame: 95x62x25cm
- Seat rest: 48x54x10cm
- Back rest: 76x54x10cm
- Leg rest: 45x54x10cm
- Arm rest: 48x15x4,5cm

Distributor:



Producer:  
VIA ACTUAL-COMÉRCIO INTERNACIONAL, LDA  
CE Castelo da Maia  
Rua Manuel Assunção Falcão, 481- Fracção H  
4475-041 Santa Maria de Avioso - Maia  
PORTUGAL

Tlf. / Phone: +351 229 829 663/4  
Fax: +351 229 829 665  
E-mail: comercial@actualway.pt  
www.actualway.pt